



Golden Silicon-Nonconductive-Non-Contact Mode- AFM Probes-c

GS-NNCM-AP-c

Our Golden silicon AFM probes have extremely sharp tips for high-resolution imaging and come standard with Au reflective coating. The standard chip sizes make these AFM probes compatible with most AFM manufacturers.

SPECIFICATIONS

short_desc	Golden Silicon AFM Probes for Nonconductive Non-Contact Mode
Categories	AFM Probes, Golden
Tags	: Contact, General Topography, Liquid Scanning, Mechanical Properties/Force Curves, Nanoindentation and Lithography, Non-Contact/Tapping, Phase Imaging Mode
Quantity	15, 50, HW, W
Cantilever Length	95
Cantilever Width	30
Cantilever Thickness	2
Resonant Frequency [kHz]	140, 240, 390
Force Constant [N/m]	3.1, 11.8, 37.6
Coating	Au Reflective
Cantilever Shape	Rectangular
Cantilever	Single
Tip	Tetrahedral (Standard)
Tip Material	Silicon

Shop online at afm-nanoprobes.cd-bioparticles.net or contact us at:

Email: info@cd-bioparticles.com | Tel: 1-631-624-4882 | Fax: 1-631-938-8221